


Simple method for the fractionation of natural compounds from plant extracts by low-resolution silica gel chromatography

Equipment and reagents required

- Vacuum concentrating evaporator
- Small clamp stand
- Clean, small bore glass pasteur pipettes (e.g. Fisher 230 mm length Pasteur pipettes, 11566963)
- Cotton buds
- 5 mg freeze-dried powder of aqueous extract, or resin of fully evaporated di-chloromethane extract of natural product of interest
- Silica 60 gel (e.g. Sigma 236799-100G)
- Acid washed sand (e.g. Sigma 18649-1KG)
- 1.5 ml microtubes or glass vials for fraction collection

Method

- 1) Suspend a clean glass pipette vertically in a clamp-stand (see right)
- 2) Insert a small piece of cotton bud from the top opening of the pipette, and tamp gently into the bottom of the column using a long piece of wire
- 3) Using a folded paper funnel or similar, layer ~3-5 mm of acid-treated sand on top of the cotton
- 4) Using a folded paper funnel or similar, layer silica 60 gel on top of the acid sand to a height reaching approximately two thirds of the way up the column
- 5) Wet the silica gel in the column by pipetting into the top opening 1 ml of the first solvent to be used for separation (see below), and discard the eluate
- 6) Prepare a slurry of ~5 mg dried extract powder (or resin from complete evaporation of a di-chloromethane extract) by mixing with a small quantity of silica gel and the first solvent
- 7) Layer this slurry on top of the silica gel already in the column
- 8) Layer a small amount of acid-treated sand on top of the slurry
- 9) Prepare sufficient solvents (see below) to perform a stepwise polarity gradient elution
- 10) Pour 2 ml of the first solvent down the column and discard the eluate
- 11) Pour 1.4 ml of each solvent down the column, in order of increasing polarity, twice for each solvent
- 12) Allow each solvent to drip through by gravity alone and collect 1.4 ml of each eluate in numbered 1.5 ml microtubes
- 13) Discard the column in a glass waste bin
- 14) Dry the eluates in each microtube using a vacuum evaporator (continue the drying process until no trace of solvent remains)
- 15) Resuspend each dried extract in 0.2 ml DMSO with vortexing
- 16) Re-assay each fraction at a maximum final concentration of 1% DMSO



Suggested solvents,* in order of increasing polarity:

- 1) 10% ethyl acetate in hexane (or octane)
- 2) 50% ethyl acetate in hexane (or octane)
- 3) 100% ethyl acetate
- 4) 10% methanol in ethyl acetate
- 5) 50% methanol in ethyl acetate
- 6) 100% methanol
- 7) 10% water in methanol
- 8) 50% water in methanol

* Note that the list of solvents suggested in this method is not definitive, and alternatives with similar step-changes in polarity would be equally suitable. Likewise, the solvents used for separation can be simplified if desired (e.g. to just 4 intermediate polarity solvents).

If desired, larger quantities of extract can be separated using larger diameter columns (e.g. 50 mg extract can be separated on a 1-2 cm diameter glass column), with volumes scaled up accordingly.