

## Method Sheet 31

# Measuring enhancement of sensitivity to existing antibiotics

### Overview

This method sheet explains how to test whether your hit extract is capable of enhancing the capacity of existing antibiotics to kill the bacterial strain you have been studying. This is a potentially useful approach for therapy, as there is much effort underway currently to find ways of defeating bacterial mechanisms of antibiotic resistance, which may in turn restore the clinical effectiveness of existing antibiotics.

### Reagents

- One or more hit natural extracts at 10 mg/ml in DMSO
- An established antibiotic compound (e.g. ampicillin, tetracycline) at 10 mg/ml

### Equipment

- Multichannel pipette (8- or 12-channel) capable of dispensing 1 - 10  $\mu$ l
- Sterile pipette tips (compatible with each pipette)
- Static (non-shaking) 37°C incubator
- Waste container for used tips

### Preparation of dilutions

- 1) Fully defrost one of the hit natural extracts and one aliquot of an antibiotic of interest (which should be at 10 mg/ml).
- 2) Prepare a 96-well plate containing bacterial cell cultures at 80  $\mu$ l per well (follow Method Sheet 03, this should be done early in the morning on the same day as challenge).
- 3) Label 10 sterile microtubes with the antibiotic concentrations 0, 0.5, 1, 2, 4, 8, 16, 32, 64 and 128.
- 4) Add 50  $\mu$ l of LB medium into each of these tubes, except the tube labelled 128, into which you should pipette 87  $\mu$ l of LB medium.
- 5) Supplement the tube labelled 128 with 12.8  $\mu$ l of antibiotic and mix by pipetting gently.
- 6) Prepare doubling dilutions of the antibiotic by transferring 50  $\mu$ l of suspension from the 128 tube into the 64 tube, mix by pipetting gently, then repeat this process by transferring 50  $\mu$ l of the resulting suspension into the next tube in the row.
- 7) Continue preparing the series of dilutions in this way until the tube labelled 0.5, but do not add anything to the tube labelled 0, this will be the zero antibiotic control.
- 8) Now label a fresh microtube with 'DMSO' and another with 'Extract'.
- 9) Pipette 25  $\mu$ l of DMSO into 225  $\mu$ l of LB medium in the 'DMSO' tube.
- 10) Pipette 25  $\mu$ l of hit extract into 225  $\mu$ l of LB medium in the 'Extract' tube.

## Challenging the bacterial cultures

- 1) Supplement all of the 20 wells shown in the plate map below with 10  $\mu$ l of the 'DMSO' control suspension (the first two rows).
- 2) Then supplement all of the 20 wells shown in the plate map below with 10  $\mu$ l of the 'Extract' suspension (the second two rows).
- 3) Supplement each column of 4 wells in the treated section with 10  $\mu$ l per well of the appropriate antibiotic dilution, as shown by the numbers in the plate map below:

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B		DMSO 0	DMSO 0.5	DMSO 1	DMSO 2	DMSO 4	DMSO 8	DMSO 16	DMSO 32	DMSO 64	DMSO 128	
C		DMSO 0	DMSO 0.5	DMSO 1	DMSO 2	DMSO 4	DMSO 8	DMSO 16	DMSO 32	DMSO 64	DMSO 128	
D		Extract 0	Extract 0.5	Extract 1	Extract 2	Extract 4	Extract 8	Extract 16	Extract 32	Extract 64	Extract 128	
E		Extract 0	Extract 0.5	Extract 1	Extract 2	Extract 4	Extract 8	Extract 16	Extract 32	Extract 64	Extract 128	
F												
G		Values indicate final concentrations of antibiotic in $\mu$ g/ml										
H												

- 4) Note that we are avoiding use of the outer wells of the plate to minimise errors arising from the 'edge effect'.
- 5) The other wells on the plate receive no supplement (i.e. they remain only 80  $\mu$ l culture).
- 6) Move the culture plate backwards and forwards, then side to side, gently several times to mix the suspension into each culture.
- 7) Avoid jarring movements that may splash liquid between wells and cause cross-contamination.
- 8) Measure the absorbance of the plate straight after challenge (to obtain the baseline, time = 0 h readings) using a microplate reader (follow Method Sheet 05).
- 9) Place the culture plate in a static 37°C incubator to allow growth to begin (does not need to be a shaking incubator).
- 10) Measure the growth once more on day 2, at 24 hours after the first measurement (i.e. at the same time as the first measurement was taken the day before).
- 11) The plate can now be discarded.
- 12) Repeat this experiment a further three times in separate plates to obtain four experimental replicates.
- 13) Proceed to analyse the data as shown in Method Sheet 32.

## Notes

- You can choose to test any antibiotic you have available in the laboratory, but good options to try in which enhancement can often be observed are tetracycline and chloramphenicol.
- This protocol suggests measuring growth at a single time point (24 h), but it will be possible to measure the absorbance of the plate at multiple timepoints to generate more data from the same experiment.
- This protocol yields a final concentration of DMSO of 1% vol/vol and a final concentration of extract at 100 µg/ml.

**Disclaimer:** These method sheets and other resources are provided for educational purposes only. The user's University Supervisor remains the Principal Investigator and the sole party responsible for the safe conduct, risk assessment, and ethical oversight of all laboratory work. Caithness Biotechnologies Ltd. accepts no liability for any injury, loss, or damage resulting from the application of the advice or protocols provided herein. Copyright © 2026, Caithness Biotechnologies Ltd. All Rights Reserved.